



# COLD SNAP™



**STEP 1**  
Soak  
with water



**STEP 2**  
Wring out  
excess water



**STEP 3**  
Snap & Shake  
to activate



ABSORBENT



LONG-LASTING



REUSABLE



WASHABLE



ANTI-MICROBIAL



CT100

CT200

CT300

CT400

CT500

## COLD SNAP™ Cooling Towels

- 33.5"L x 13"W
- Absorbent PVC material accelerates cooling process
- Stays cool for hours
- Reusable: Rinse, wring, fold & store in tube
- Machine or hand wash, air dry
- Anti-microbial treated to prevent mildew & odor



CB100

CB200

CB300

CB400

CB500

## COLD SNAP™ Cooling Bands

- 29.4"L x 4"W
- Absorbent PVC material accelerates cooling process
- Stays cool for hours
- Reusable: Rinse, wring, fold & store in tube
- Machine or hand wash, air dry
- Anti-microbial treated to prevent mildew & odor



## COLD SNAP™ Cooling Head Bands | #CHB100

- 34"L x 2"W, 2-Pack
- Absorbent PVC material accelerates cooling process
- Stays cool for hours
- Reusable: Rinse, wring, fold & store in tube
- Machine or hand wash, air dry
- Anti-microbial treated to prevent mildew & odor